

Antibodies Product Datasheet

Product Name:	MYL3 Monoclonal Antibody
Catalog Number:	N0458
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of MYL3 expressed in E. Coli.
Concentration:	
Gene Name:	MYL3
Protein Name:	Myosin light chain 3
Other Name:	MYL3; Myosin light chain 3; Cardiac myosin light chain 1; CMLC1; Myosin light chain 1; slow-twitch muscle B/ventricular isoform; MLC1SB; Ventricular/slow twitch myosin alkali light chain
Human Swiss Prot No.:	P08590
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	MYL3 Monoclonal Antibody detects endogenous levels of MYL3 protein.
Storage and Stability:	-20°C/1 year