

Antibodies Product Datasheet

Product Name:	mGluR7 Antibody
Catalog Number:	G367
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRM7. AA range:351-400
Concentration:	1 mg/ml
Gene Name:	GRM7
Protein Name:	Metabotropic glutamate receptor 7
Other Name:	GRM7; GPRC1G; MGLUR7; Metabotropic glutamate receptor 7; mGluR7
Human Swiss Prot No.:	Q14831
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	mGluR7 Polyclonal Antibody detects endogenous levels of mGluR7 protein.
Storage and Stability:	-20°C/1 year