

Antibodies Product Datasheet

Product Name:	YY1 (Acetyl-Lys351) Antibody
Catalog Number:	D12222
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human YY1 (Acetyl Lys351)
Concentration:	1 mg/ml
Gene Name:	YY1 INO80S
Protein Name:	YY1 (Acetyl Lys351)
Other Name:	Transcriptional repressor protein YY1 (Delta transcription factor) (INO80 complex subunit S) (NF-E1) (Yin and yang 1) (YY-1)
Human Swiss Prot No.:	P25490
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse YY1 (Acetyl Lys351)
Storage and Stability:	-20°C/1 year