

## Antibodies Product Datasheet

<b>Product Name:</b>	TIS11B Antibody
<b>Catalog Number:</b>	B1184
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human TISB. AA range:58-107
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ZFP36L1
<b>Protein Name:</b>	Zinc finger protein 36 C3H1 type-like 1
<b>Other Name:</b>	ZFP36L1; BERG36; BRF1; ERF1; RNF162B; TIS11B; Zinc finger protein 36; C3H1 type-like 1; Butyrate response factor 1; EGF-response factor 1; ERF-1; Protein TIS11B
<b>Human Swiss Prot No.:</b>	Q07352
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	TIS11B Polyclonal Antibody detects endogenous levels of TIS11B protein.
<b>Storage and Stability:</b>	-20°C/1 year