

Antibodies Product Datasheet

Product Name:	CASC5/KNL1 Antibody
Catalog Number:	C51078
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CASC5/KNL1
Concentration:	1 mg/ml
Gene Name:	CASC5 KIAA1570 KNL1
Protein Name:	Protein CASC5 (ALL1-fused gene from chromosome 15q14 protein) (AF15q14) (Bub-linking kinetochore protein) (Blinkin) (Cancer susceptibility candidate gene 5 protein) (Cancer/testis antigen 29) (CT29) (Kinetochore-null protein 1) (Protein D40/AF15q14)
Other Name:	
Human Swiss Prot No.:	Q8NG31
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CASC5/KNL1 at Human
Storage and Stability:	Store at -20° C/1 year