

Antibodies Product Datasheet

Product Name:	Cdc16 Antibody
Catalog Number:	C12166
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APC6. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	CDC16
Protein Name:	Cell division cycle protein 16 homolog
Other Name:	CDC16; ANAPC6; Cell division cycle protein 16 homolog; Anaphase-promoting complex subunit 6; APC6; CDC16 homolog; CDC16Hs; Cyclosome subunit 6
Human Swiss Prot No.:	Q13042
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc16 Polyclonal Antibody detects endogenous levels of Cdc16 protein.
Storage and Stability:	-20°C/1 year