

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Rabbit Anti Sheep IgG (H&L) - Alexa Fluor 680
<b>Catalog Number:</b>	SA1028
<b>Reactivity:</b>	Sheep
<b>Conjugate:</b>	Alexa Fluor 680
<b>Applications:</b>	IF/ICC,FCM,WB
<b>Dilutions:</b>	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Concentration:</b>	1mg/mL
<b>Other Name:</b>	Rabbit Anti-Sheep IgG (H&L) - Alexa Fluor 680
<b>Formulation:</b>	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
<b>Storage and Stability:</b>	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing