

## Antibodies Product Datasheet

<b>Product Name:</b>	Dipeptidyl-peptidase 1 (heavy chain, Cleaved-Leu231) Antibody
<b>Catalog Number:</b>	L0460
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Modification:</b>	Cleaved Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human Dipeptidyl-peptidase 1 (heavy chain, Cleaved-Leu231)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CTSC CPPI
<b>Protein Name:</b>	Dipeptidyl-peptidase 1 (heavy chain, Cleaved-Leu231)
<b>Other Name:</b>	Dipeptidyl peptidase 1 (EC 3.4.14.1) (Cathepsin C) (Cathepsin J) (Dipeptidyl peptidase I) (DPP-I) (DPPI) (Dipeptidyl transferase) [Cleaved into: Dipeptidyl peptidase 1 exclusion domain chain (Dipeptidyl peptidase I exclusion domain chain); Dipeptidyl pept
<b>Human Swiss Prot No.:</b>	P53634
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human Dipeptidyl-peptidase 2 (heavy chain, Cleaved-Leu231, protein was cleaved amino acid sequence between 230-231 )
<b>Storage and Stability:</b>	-20°C/1 year