

Antibodies Product Datasheet

Product Name:	Anti-Oct-2 [ABT-OCT2]
Catalog Number:	V0105
Clone ID:	ABT-OCT2
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	P09086
Gene Name:	POU2F2
Other Name:	POU domain class 2 transcription factor 2; Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2; Octamer-binding protein 2; Oct-2; Octamer-binding transcription factor 2; OTF-2 ; POU2F2 ; OCT2; OTF2
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Oct-2 [ABT-OCT2] detects endogenous levels of total Oct-2 protein.
Storage and Stability:	-20°C/1 year