

Antibodies Product Datasheet

Product Name:	Cytokeratin 10 Antibody
Catalog Number:	C0249
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Keratin 10. AA range:136-185
Concentration:	1 mg/ml
Gene Name:	KRT10
Protein Name:	Keratin type I cytoskeletal 10
Other Name:	KRT10; KPP; Keratin; type I cytoskeletal 10; Cytokeratin-10; CK-10; Keratin-10; K10
Human Swiss Prot No.:	P13645
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cytokeratin 10 Polyclonal Antibody detects endogenous levels of Cytokeratin 10 protein.
Storage and Stability:	-20°C/1 year