

Antibodies Product Datasheet

Product Name:	RHBDD1 Antibody
Catalog Number:	C48661
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RHBDD1
Concentration:	1 mg/ml
Gene Name:	RHBDD1 RHBDL4 HSD-50 HSD50
Protein Name:	Rhomboid-related protein 4 (RRP4) (EC 3.4.21.105) (Rhomboid domain-containing protein 1) (Rhomboid-like protein 4)
Other Name:	
Human Swiss Prot No.:	Q8TEB9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RHBDD1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year