

Antibodies Product Datasheet

Product Name:	NGFR p75 Antibody
Catalog Number:	R12-2265
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NGFR p75. AA range:230-279
Concentration:	1 mg/ml
Gene Name:	NGFR
Protein Name:	Tumor necrosis factor receptor superfamily member 16
Other Name:	NGFR; TNFRSF16; Tumor necrosis factor receptor superfamily member 16; Gp80-LNGFR; Low affinity neurotrophin receptor p75NTR; Low-affinity nerve growth factor receptor; NGF receptor; p75 ICD; CD antigen CD271
Human Swiss Prot No.:	P08138
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NGFR p75 Polyclonal Antibody detects endogenous levels of NGFR p75 protein.
Storage and Stability:	-20°C/1 year