
Antibodies Product Datasheet

Product Name:	CCDC102A Antibody
Catalog Number:	C15185
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCDC102A. AA range:341-390
Concentration:	1 mg/ml
Gene Name:	CCDC102A
Protein Name:	Coiled-coil domain-containing protein 102A
Other Name:	CCDC102A; Coiled-coil domain-containing protein 102A
Human Swiss Prot No.:	Q96A19
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CCDC102A Polyclonal Antibody detects endogenous levels of CCDC102A protein.
Storage and Stability:	-20°C/1 year