

Antibodies Product Datasheet

Product Name:	FMO2 Antibody
Catalog Number:	C50494
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FMO2
Concentration:	1 mg/ml
Gene Name:	FMO2
Protein Name:	Dimethylaniline monooxygenase [N-oxide-forming] 2 (EC 1.14.13.8) (Dimethylaniline oxidase 2) (FMO 1B1) (Pulmonary flavin-containing monooxygenase 2) (FMO 2)
Other Name:	
Human Swiss Prot No.:	Q99518
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FMO2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year