

Antibodies Product Datasheet

Product Name:	Cryopyrin Antibody
Catalog Number:	C30251
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human NLRP3. AA range:511-560
Concentration:	1 mg/ml
Gene Name:	NLRP3
Protein Name:	NACHT LRR and PYD domains-containing protein 3
Other Name:	NLRP3; C1orf7; CIAS1; NALP3; PYPAF1; NACHT, LRR and PYD domains-containing protein 3; Angiotensin/vasopressin receptor AII/AVP-like; Caterpillar protein 1.1CLR1.1; Cold autoinflammatory syndrome 1 protein; Cryopyrin; PYRIN-containing APAF1-like protein 1
Human Swiss Prot No.:	Q96P20
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cryopyrin Polyclonal Antibody detects endogenous levels of Cryopyrin protein.
Storage and Stability:	-20°C/1 year