

Antibodies Product Datasheet

Product Name:	CD161 Antibody
Catalog Number:	C30502
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human KLRB1. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	KLRB1
Protein Name:	Killer cell lectin-like receptor subfamily B member 1
Other Name:	KLRB1; CLEC5B; NKRP1A; Killer cell lectin-like receptor subfamily B member 1; C-type lectin domain family 5 member B; HNKR-P1a; NKRP1A; Natural killer cell surface protein P1A; CD161
Human Swiss Prot No.:	Q12918
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD161 Polyclonal Antibody detects endogenous levels of CD161 protein.
Storage and Stability:	-20°C/1 year