
Antibodies Product Datasheet

Product Name:	Cyclin A Antibody
Catalog Number:	C0165
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cyclin A. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	CCNA1/CCNA2
Protein Name:	Cyclin-A1/2
Other Name:	CCNA1; Cyclin-A1; CCNA2; CCN1; CCNA; Cyclin-A2; Cyclin-A
Human Swiss Prot No.:	P78396/P20248
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclin A Polyclonal Antibody detects endogenous levels of Cyclin A protein.
Storage and Stability:	-20°C/1 year