

Antibodies Product Datasheet

Product Name:	RIAM Antibody
Catalog Number:	C14452
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from RIAM . at AA range: 430-510
Concentration:	1 mg/ml
Gene Name:	APBB1IP
Protein Name:	Amyloid beta A4 precursor protein-binding family B member 1-interacting protein
Other Name:	APBB1IP; PREL1; RARP1; RIAM; Amyloid beta A4 precursor protein-binding family B member 1-interacting protein; APBB1-interacting protein 1; Proline-rich EVH1 ligand 1; PREL-1; Proline-rich protein 73; Rap1-GTP-interacting adapter molecule; R
Human Swiss Prot No.:	Q7Z5R6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RIAM Polyclonal Antibody detects endogenous levels of RIAM protein.
Storage and Stability:	-20°C/1 year