

Antibodies Product Datasheet

Product Name:	Bek Antibody
Catalog Number:	C0191
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FGFR2. AA range:471-520
Concentration:	1 mg/ml
Gene Name:	FGFR2
Protein Name:	Fibroblast growth factor receptor 2
Other Name:	FGFR2; BEK; KGFR; KSAM; Fibroblast growth factor receptor 2; FGFR-2; K-sam; KGFR; Keratinocyte growth factor receptor; CD antigen CD332
Human Swiss Prot No.:	P21802
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Bek Polyclonal Antibody detects endogenous levels of Bek protein.
Storage and Stability:	-20°C/1 year