

Antibodies Product Datasheet

Product Name:	ABCG1 Antibody
Catalog Number:	C14632
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ABCG1 . at AA range: 560-640
Concentration:	1 mg/ml
Gene Name:	ABCG1
Protein Name:	ATP-binding cassette sub-family G member 1
Other Name:	ABCG1; ABC8; WHT1; ATP-binding cassette sub-family G member 1; ATP-binding cassette transporter 8; White protein homolog
Human Swiss Prot No.:	P45844
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABCG1 Polyclonal Antibody detects endogenous levels of ABCG1 protein.
Storage and Stability:	-20°C/1 year