
Antibodies Product Datasheet

Product Name:	FoxO4 Antibody
Catalog Number:	B8037
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AFX1. AA range:228-277
Concentration:	1 mg/ml
Gene Name:	FOXO4
Protein Name:	Forkhead box protein O4
Other Name:	FOXO4; AFX; AFX1; MLLT7; Forkhead box protein O4; Fork head domain transcription factor AFX1
Human Swiss Prot No.:	P98177
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FoxO4 Polyclonal Antibody detects endogenous levels of FoxO4 protein.
Storage and Stability:	-20°C/1 year