

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Histone H4 (Mono Methyl Arg3) Antibody |
| Catalog Number: | D12834 |
| Reactivity: | Human |
| Applications: | WB |
| Dilutions: | WB: 1:1000-2000 |
| Modification: | Methyl Specific |
| Source: | Polyclonal Mouse |
| Purification: | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthetic Peptide of Histone H4 (Mono Methyl Arg3) |
| Concentration: | |
| Gene Name: | HIST1H4A |
| Protein Name: | Histone H4 |
| Other Name: | HIST1H4A; H4/A; H4FA; HIST1H4B; H4/I; H4FI; HIST1H4C; H4/G; H4FG; HIST1H4D; H4/B; H4FB; HIST1H4E; H4/J; H4FJ; HIST1H4F; H4/C; H4FC; HIST1H4H; H4/H; H4FH; HIST1H4I; H4/M; H4FM; HIST1H4J; H4/E; H4FE; HIST1H4K; H4/D; H4FD; HIST1H4L; H4/K; H4FK; HIST2H4A; H4/ |
| Human Swiss Prot No.: | P62805 |
| Formulation: | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Specificity: | The antibody detects endogenous Histone H3 (Mono Methyl Arg3) protein. |
| Storage and Stability: | -20°C/1 year |