

---

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | BIOTIN*Polyclonal Goat Anti Human IgG Fab   |
| <b>Catalog Number:</b>        | SA0306  |
| <b>Reactivity:</b>            | HumanIgGFab   |
| <b>Conjugate:</b>             | Biotin  |
| <b>Applications:</b>          | WB,IHC-p,IF/ICC,ELISA   |
| <b>Dilutions:</b>             | WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000                          |
| <b>Modification:</b>          | Not Modification Specific   |
| <b>Source:</b>                | Goat  |
| <b>Purification:</b>          | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| <b>Concentration:</b>         | 4-5mg/mL  |
| <b>Other Name:</b>            | BIOTIN*Polyclonal Goat Anti-Human IgG Fab   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                         |
| <b>Storage and Stability:</b> | -20°C/1 year  |