

Antibodies Product Datasheet

Product Name:	Claudin 3 (PT1687) mouse mAb
Catalog Number:	N1829
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	WB 500-2000 IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Claudin 3
Concentration:	
Gene Name:	CLDN3 C7orf1 CPETR2
Protein Name:	Claudin-3 (Clostridium perfringens enterotoxin receptor 2) (CPE-R 2) (CPE-receptor 2) (Rat ventral prostate.1 protein homolog) (hRVP1)
Other Name:	
Human Swiss Prot No.:	O15551
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Claudin 3. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year