
Antibodies Product Datasheet

Product Name:	GATA4 Antibody
Catalog Number:	R12-2799
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	GATA4 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GATA4.
Concentration:	1 mg/ml
Gene Name:	GATA4
Other Name:	DNA binding protein GATA-GT2; GAT4; GATA binding factor-4; GATA-4; Transcription factor GATA-4
Human Swiss Prot No.:	P43694
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	GATA4 Antibody antibody detects endogenous levels of GATA4.
Storage and Stability:	Stable at -20°C for at least 1 year.