

## Antibodies Product Datasheet

<b>Product Name:</b>	CD3EAP Antibody
<b>Catalog Number:</b>	R12-2588
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	CD3EAP Antibody
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Rabbit Polyclonal
<b>Purification:</b>	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human CD3EAP.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CD3EAP
<b>Other Name:</b>	DNA-directed RNA polymerase I subunit RPA34; RNA polymerase I-associated factor PAF49; Anti-sense to ERCC-1 protein; ASE-1; CD3-epsilon-associated protein; CD3E-associated protein; CAST; A34.5; ASE1; CAST; PAF49; CD3EAP
<b>Human Swiss Prot No.:</b>	O15446
<b>Formulation:</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Specificity:</b>	CD3EAP Antibody antibody detects endogenous levels of CD3EAP.
<b>Storage and Stability:</b>	Stable at -20°C for at least 1 year.