

Antibodies Product Datasheet

Product Name:	HP-1 α mouse Monoclonal Antibody(5E3)
Catalog Number:	N1687
Reactivity:	Human,Mouse,Rat
Applications:	IF/ICC,WB,IHC-p
Dilutions:	IF: 1:50-200 WB 1:500-2000,IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of HP-1 α
Concentration:	1 mg/ml
Gene Name:	CBX5 HP1A
Protein Name:	Chromobox protein homolog 5 (Antigen p25) (Heterochromatin protein 1 homolog alpha) (HP1 alpha)
Other Name:	Chromobox protein homolog 5 (Antigen p25) (Heterochromatin protein 1 homolog alpha) (HP1 alpha)
Human Swiss Prot No.:	P45973
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous HP-1 α protein
Storage and Stability:	-20°C/1 year