

Antibodies Product Datasheet

Product Name:	TARSL2 Antibody
Catalog Number:	C50825
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TARSL2
Concentration:	1 mg/ml
Gene Name:	TARSL2
Protein Name:	Probable threonine--tRNA ligase 2, cytoplasmic (EC 6.1.1.3) (Threonyl-tRNA synthetase) (ThrRS) (Threonyl-tRNA synthetase-like protein 2)
Other Name:	
Human Swiss Prot No.:	A2RTX5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TARSL2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year