
Antibodies Product Datasheet

Product Name:	XRCC4 Monoclonal Antibody(5C10)
Catalog Number:	N1351
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCIP
Dilutions:	WB: 1:2000 IP:1:200 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of XRCC4
Concentration:	
Gene Name:	XRCC4
Protein Name:	DNA repair protein XRCC4
Other Name:	XRCC4; DNA repair protein XRCC4; X-ray repair cross-complementing protein 4
Human Swiss Prot No.:	Q13426
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous XRCC4 proteins.
Storage and Stability:	-20°C/1 year