

Antibodies Product Datasheet

Product Name:	TMBIM1 Antibody
Catalog Number:	R12-3629
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	TMBIM1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TMBIM1.
Concentration:	1 mg/ml
Gene Name:	TMBIM1
Other Name:	Transmembrane BAX inhibitor motif-containing protein 1; Protein RECS1 homolog; TMBIM1; RECS1; PP1201; PSEC0158
Human Swiss Prot No.:	Q969X1
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TMBIM1 Antibodyantibody detects endogenous levels of TMBIM1.
Storage and Stability:	Stable at -20°C for at least 1 year.