

Antibodies Product Datasheet

Product Name:	COX IV Monoclonal Antibody(6C8)
Catalog Number:	N1297
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:1000-3000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Concentration:	
Gene Name:	COX4I1
Protein Name:	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Other Name:	COX4I1; COX4; Cytochrome c oxidase subunit 4 isoform 1, mitochondrial; Cytochrome c oxidase polypeptide IV; Cytochrome c oxidase subunit IV isoform 1; COX IV-1
Human Swiss Prot No.:	P13073
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous COX IV protein.
Storage and Stability:	-20°C/1 year