

Antibodies Product Datasheet

Product Name:	B-Myb (Phospho-Ser577) Antibody
Catalog Number:	A0823
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human B-Myb around the phosphorylation site of Ser577/581. AA range:551-600
Concentration:	1 mg/ml
Gene Name:	MYBL2
Protein Name:	Myb-related protein B
Other Name:	MYBL2; BMYB; Myb-related protein B; B-Myb; Myb-like protein 2
Human Swiss Prot No.:	P10244
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-B-Myb (S577) Polyclonal Antibody detects endogenous levels of B-Myb protein only when phosphorylated at S577.
Storage and Stability:	-20°C/1 year