

Antibodies Product Datasheet

Product Name:	NR2F1 Antibody
Catalog Number:	C48937
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NR2F1
Concentration:	1 mg/ml
Gene Name:	NR2F1 EAR3 ERBAL3 TFCOUP1
Protein Name:	COUP transcription factor 1 (COUP-TF1) (COUP transcription factor I) (COUP-TF I) (Nuclear receptor subfamily 2 group F member 1) (V-erbA-related protein 3) (EAR-3)
Other Name:	
Human Swiss Prot No.:	P10589
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NR2F1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year