

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | INMT Antibody |
| Catalog Number: | C50259 |
| Reactivity: | Human, Mouse |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit, IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human INMT |
| Concentration: | 1 mg/ml |
| Gene Name: | INMT |
| Protein Name: | Indolethylamine N-methyltransferase (Indolamine N-methyltransferase) (EC 2.1.1.49) (EC 2.1.1.96) (Aromatic alkylamine N-methyltransferase) (Amine N-methyltransferase) (Arylamine N-methyltransferase) (Thioether S-methyltransferase) (TEMT) |
| Other Name: | |
| Human Swiss Prot No.: | O95050 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of INMT at Human, Mouse |
| Storage and Stability: | Store at -20° C/1 year |