

Antibodies Product Datasheet

Product Name:	LACTB Antibody
Catalog Number:	C48405
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LACTB
Concentration:	1 mg/ml
Gene Name:	LACTB MRPL56 UNQ843/PRO1781
Protein Name:	Serine beta-lactamase-like protein LACTB, mitochondrial (EC 3.4.-.-)
Other Name:	
Human Swiss Prot No.:	P83111
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LACTB at Human, Mouse
Storage and Stability:	Store at -20° C/1 year