

Antibodies Product Datasheet

Product Name:	Neuregulin-1 SMDF Antibody
Catalog Number:	C0229
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NRG1 isoform-10. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	NRG1
Protein Name:	Pro-neuregulin-1 membrane-bound isoform
Other Name:	Pro-neuregulin-1 membrane-bound isoform; Pro-NRG1; NRG1; GGF; HGL; HRGA; NDF; SMDF
Human Swiss Prot No.:	Q02297-10
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Neuregulin-1 SMDF Polyclonal Antibody detects endogenous levels of Neuregulin-1 SMDF protein.
Storage and Stability:	-20°C/1 year