

Antibodies Product Datasheet

Product Name:	MIF Antibody
Catalog Number:	R12-2238
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MIF. AA range:25-74
Concentration:	1 mg/ml
Gene Name:	MIF
Protein Name:	Macrophage migration inhibitory factor
Other Name:	MIF; GLIF; MMIF; Macrophage migration inhibitory factor; MIF; Glycosylation-inhibiting factor; GIF; L-dopachrome isomerase; L-dopachrome tautomerase; Phenylpyruvate tautomerase
Human Swiss Prot No.:	P14174
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MIF Polyclonal Antibody detects endogenous levels of MIF protein.
Storage and Stability:	-20°C/1 year