

Antibodies Product Datasheet

Product Name:	ZC3H15 Antibody
Catalog Number:	C50681
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ZC3H15
Concentration:	1 mg/ml
Gene Name:	ZC3H15 DFRP1 LEREPO4 HSPC303 HT010 MSTP012 PP730
Protein Name:	Zinc finger CCCH domain-containing protein 15 (DRG family-regulatory protein 1) (Likely ortholog of mouse immediate early response erythropoietin 4)
Other Name:	
Human Swiss Prot No.:	Q8WU90
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ZC3H15 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year