
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | STK39 Antibody |
| Catalog Number: | R12-3572 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | STK39 Antibody |
| Modification: | Unmodified/Total |
| Source: | Rabbit Polyclonal |
| Purification: | The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human STK39. |
| Concentration: | 1 mg/ml |
| Gene Name: | STK39 |
| Other Name: | Ste-20-related kinase; Serine/threonine-protein kinase 39; DCHT |
| Human Swiss Prot No.: | Q9UEW8 |
| Formulation: | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Specificity: | STK39 Antibody antibody detects endogenous levels of STK39. |
| Storage and Stability: | Stable at -20°C for at least 1 year. |