

## Antibodies Product Datasheet

<b>Product Name:</b>	PAR15 Antibody
<b>Catalog Number:</b>	C45417
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human protein . at AA range: 170-250
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PARP15 BAL3
<b>Protein Name:</b>	Poly [ADP-ribose] polymerase 15 (PARP-15) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 7) (ARTD7) (B-aggressive lymphoma protein 3)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q460N3
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	PAR15 Polyclonal Antibody detects endogenous levels of protein.
<b>Storage and Stability:</b>	-20°C/1 year