

Antibodies Product Datasheet

Product Name:	CD71 Antibody
Catalog Number:	R12-2082
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	CD71 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD71.
Concentration:	1 mg/ml
Gene Name:	TFRC
Other Name:	Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD71; Transferrin receptor protein 1, serum form; sTfR; GTFRC
Human Swiss Prot No.:	P02786
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CD71 antibody detects endogenous levels of CD71.
Storage and Stability:	Stable at -20°C for at least 1 year.