

Antibodies Product Datasheet

Product Name:	MYL6 Antibody
Catalog Number:	C20889
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human MYL6.
Concentration:	1 mg/ml
Gene Name:	MYL6
Protein Name:	Myosin light polypeptide 6
Other Name:	MYL6; Myosin light polypeptide 6; 17 kDa myosin light chain; LC17; Myosin light chain 3; MLC-3; Myosin light chain alkali 3; Myosin light chain A3; Smooth muscle and nonmuscle myosin light chain alkali 6
Human Swiss Prot No.:	P60660
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MYL6 Polyclonal Antibody detects endogenous levels of MYL6 protein.
Storage and Stability:	-20°C/1 year