

## Antibodies Product Datasheet

<b>Product Name:</b>	HSP40 Mouse mAb(Mix-mA)
<b>Catalog Number:</b>	N1460
<b>Reactivity:</b>	Human, Mouse,Rat
<b>Applications:</b>	IHC-p,WB
<b>Dilutions:</b>	IHC-p1:50-200 ,WB 1:1000-2000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human HSP40
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	DNAJB1 DNAJ1 HDJ1 HSPF1
<b>Protein Name:</b>	HSP40
<b>Other Name:</b>	DnaJ homolog subfamily B member 1 (DnaJ protein homolog 1) (Heat shock 40 kDa protein 1) (HSP40) (Heat shock protein 40) (Human DnaJ protein 1) (hDj-1)
<b>Human Swiss Prot No.:</b>	P25685
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.33% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of HSP40 at Human, Mouse,Rat
<b>Storage and Stability:</b>	-20°C/1 year