
Antibodies Product Datasheet

Product Name:	ApoM Monoclonal Antibody
Catalog Number:	N0038
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ApoM expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	APOM
Protein Name:	Apolipoprotein M
Other Name:	APOM; G3A; NG20; HSPC336; Apolipoprotein M; Apo-M; ApoM; Protein G3a
Human Swiss Prot No.:	O95445
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	ApoM Monoclonal Antibody detects endogenous levels of ApoM protein.
Storage and Stability:	-20°C/1 year