

Antibodies Product Datasheet

Product Name:	Anti-Human Nucleophosmin (ABT358) Mouse Monoclonal Antibody
Catalog Number:	V0421
Clone ID:	ABT358
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	NPM1
Gene Name:	P06748
Other Name:	B23, MGC104254, NO38, NPM, NPM_HUMAN, NPM1, Nucleolar phosphoprotein B23, Nucleolar protein NO38, Nucleophosmin (nucleolar phosphoprotein B23 numatrin), Nucleophosmin, nucleophosmin nucleoplasmin family member 1, Nucleophosmin/nucleoplasmin family member 1, Numatrin, OTTHUMP00000161024, OTTHUMP00000161025, OTTHUMP00000223397, OTTHUMP00000223398
Isotype	IgG2b, kappa
Localization:	Liver
Positive Tissue Control:	Nuclear, Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Nucleophosmin protein.
Storage and Stability:	-20°C/1 year