

Antibodies Product Datasheet

Product Name:	GALR1 Antibody
Catalog Number:	G267
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GALR1. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	GALR1
Protein Name:	Galanin receptor type 1
Other Name:	GALR1; GALNR; GALNR1; Galanin receptor type 1; GAL1-R; GALR-1
Human Swiss Prot No.:	P47211
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GALR1 Polyclonal Antibody detects endogenous levels of GALR1 protein.
Storage and Stability:	-20°C/1 year