

Antibodies Product Datasheet

Product Name:	MLH1 mouse mAb
Catalog Number:	N0950
Reactivity:	Human
Applications:	WB
Dilutions:	wb dilution 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human MLH1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	mlh1
Protein Name:	
Other Name:	COCA 2;COCA2;DNA mismatch repair protein Mlh1;FCC 2;FCC2;hMLH1;hMLH1;HNPCC 2;HNPCC; HNPCC2;MGC5172;MLH1;MLH1;MLH1_HUMAN;MutL homolog 1 (E. coli);MutL homolog 1;MutL homolog 1 colon cancer nonpolyposis type 2;MutL homolog 1, colon cancer, nonpolyposis ty
Human Swiss Prot No.:	P40692
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MLH1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year