

Antibodies Product Datasheet

Product Name:	Nir2 Antibody
Catalog Number:	C49446
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Nir2
Concentration:	1 mg/ml
Gene Name:	PITPNM1 DRES9 NIR2 PITPNM
Protein Name:	Membrane-associated phosphatidylinositol transfer protein 1 (Drosophila retinal degeneration B homolog) (Phosphatidylinositol transfer protein, membrane-associated 1) (PITPnm 1) (Pyk2 N-terminal domain-interacting receptor 2) (NIR-2)
Other Name:	
Human Swiss Prot No.:	O00562
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Nir2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year