

Antibodies Product Datasheet

Product Name:	IL-1F8 Antibody
Catalog Number:	C30613
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human IL36B. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	IL36B
Protein Name:	Interleukin-36 beta
Other Name:	IL36B; IL1F8; IL1H2; Interleukin-36 beta; FIL1 eta; Interleukin-1 eta; IL-1 eta; Interleukin-1 family member 8; IL-1F8; Interleukin-1 homolog 2; IL-1H2
Human Swiss Prot No.:	Q9NZH7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-1F8 Polyclonal Antibody detects endogenous levels of IL-1F8 protein.
Storage and Stability:	-20°C/1 year