

Antibodies Product Datasheet

Product Name:	PTG Antibody
Catalog Number:	C18015
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPP1R3C. AA range:44-93
Concentration:	1 mg/ml
Gene Name:	PPP1R3C
Protein Name:	Protein phosphatase 1 regulatory subunit 3C
Other Name:	PPP1R3C; PPP1R5; Protein phosphatase 1 regulatory subunit 3C; Protein phosphatase 1 regulatory subunit 5; PP1 subunit R5; Protein targeting to glycogen; PTG
Human Swiss Prot No.:	Q9UQK1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PTG Polyclonal Antibody detects endogenous levels of PTG protein.
Storage and Stability:	-20°C/1 year